
Experimental setups in Montes- Hernandez Lab

General Topic: Nucleation, growth and transformation of minerals and particles under mild and hydrothermal conditions

Real-time monitoring of gas-solid and/or gas-liquidsolid reactions by Raman spectroscopy





Real-time monitoring of pH in pressured and hydrothermal systems



- Up to 100 bar
- Up to 150°C
- Duration: from minutes to weeks
- Time-pH acquisition: from 1s to hours

oH (in situ) Ex. 2. Direct steel slag carbonation using compressed carbon dioxide (initial pressure of 10 bar of CO2)






Ex. 1. Calcite formation in presence of organic additives from aqueous carbonation of portlandite at 55 bar of CO2 pressure...

Complementary hydrothermal reactors and ovens for mineral alteration studies and/or mineral synthesis

3 ovens of 240L for hydrothermal reactions up to 100°C and/or drying of solid materials


10 mini-reactors (3ml) for hydrothermal reactions up to 200°C and independent temperature for each reactor



- 2L up to 90 bar and 90°C
- 1L up to 260°C and 70 bar
- 0.3L up to 300°C and 300 bar


Mineral hydrothermal alteration in tubular reactors under water supercritical conditions





Flowthrough reactors coupled with multiionic and pH probes Room T expering pollutant sequences.

Room T experiments for pollutant sequestration and mineral precipitation...



Inverted microscope for crystallization

experiments

Microscope Olympus IX70





F-View camera

Crystallization from μ -drops ionic solutions and/or by mixing ionic solutions... at room T

Images from S. Marchaland-Le Bihan (L2 and DUT reports)



